





PRPH Antibody

Product Code	CSB-PA018774ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41219
Immunogen	Recombinant Human Peripherin protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peripherin (Neurofilament 4), PRPH, NEF4 PRPH1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PRPH

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA018774ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA018774ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.