





RAB5B Antibody

Product Code	CSB-PA019214ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61020
Immunogen	Recombinant Human Ras-related protein Rab-5B protein (1-215AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-5B, RAB5B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RAB5B

Image



Western blot
 All lanes: RAB5B antibody at 5.16µg/ml
 Lane 1: Hela whole cell lysate
 Lane 2: Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 24, 20 kDa
 Observed band size: 24 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA019214ESR1HU at dilution of 1:100