




RAB5B Antibody

Product Code	CSB-PA019214ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61020
Immunogen	Recombinant Human Ras-related protein Rab-5B protein (1-215AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-5B, RAB5B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RAB5B

Image



Western blot

All lanes: RAB5B antibody at 1.66μg/ml

Lane 1: A431 whole cell lysate

Lane 2: Hela whole cell lysate

Lane 3: Jurkat whole cell lysate

Lane 4: 293T whole cell lysate


Lane 5: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 24, 20 kDa

Observed band size: 24, 31 kDa



Immunohistochemistry of paraffin-embedded

human testis tissue using CSB-

PA019214ESR2HU at dilution of 1:100