





RFPL2 Antibody

Product Code	CSB-PA019599ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75678
Immunogen	Recombinant Human Ret finger protein-like 2 protein (99-378AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ret finger protein-like 2 (RING finger protein 79), RFPL2, RNF79
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RFPL2


Image



Western blot
 All lanes: RFPL2 antibody at 3.36µg/ml + MCF-7 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 43, 33, 36 kDa
 Observed band size: 43 kDa



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA019599ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA019599ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.