



SCAND1 Antibody

Product Code	CSB-PA020755DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57086
Immunogen	Recombinant Human SCAN domain-containing protein 1 protein (1-179AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	SCAN domain-containing protein 1, SCAND1, SDP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SCAND1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.