





SIK3 Antibody

Product Code	CSB-PA897531ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2K2
Immunogen	Recombinant Human Serine/threonine-protein kinase SIK3 protein (60-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase SIK3 (EC 2.7.11.1) (Salt-inducible kinase 3) (SIK-3) (Serine/threonine-protein kinase QSK), SIK3, KIAA0999 QSK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SIK3

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA897531ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human thymus tissue using CSB-PA897531ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.