





TBPL1 Antibody

Product Code	CSB-PA023240ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62380
Immunogen	Recombinant Human TATA box-binding protein-like protein 1 protein (1-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	TATA box-binding protein-like protein 1 (TBP-like protein 1) (21 kDa TBP-like protein) (Second TBP of unique DNA protein) (STUD) (TATA box-binding protein-related factor 2) (TBP-related factor 2) (TBP-like factor) (TBP-related protein), TBPL1, TLF TLP TLP21 TRF2 TRP
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TBPL1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA023240ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA023240ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.