




ZNF785 Antibody

Product Code	CSB-PA027041ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8K8V0
Immunogen	Recombinant Human Zinc finger protein 785 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 785, ZNF785
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF785

Image



Western blot

All lanes: ZNF785 antibody at 6.27μg/ml

Lane 1: Hela whole cell lysate


Lane 2: PC-3 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 47, 45 kDa

Observed band size: 47, 65, 30 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA027041ESR2HU at dilution of 1:100