





UGT2B10 Antibody

Product Code	CSB-PA025587ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36537
Immunogen	Recombinant Human UDP-glucuronosyltransferase 2B10 protein (20-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-glucuronosyltransferase 2B10 (UDPGT 2B10) (EC 2.4.1.17), UGT2B10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UGT2B10

Image



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA025587ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA025587ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.