




WDR45 Antibody

Product Code	CSB-PA852773ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y484
Immunogen	Recombinant Human WD repeat domain phosphoinositide-interacting protein 4 protein (1-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat domain phosphoinositide-interacting protein 4 (WIPI-4) (WD repeat-containing protein 45), WDR45, WDRX1 WDRX14 WIPI4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WDR45

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA852773ESR1HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.