





POC1A Antibody

Product Code	CSB-PA822753ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBT0
Immunogen	Recombinant Human POC1 centriolar protein homolog A protein (138-407AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	POC1 centriolar protein homolog A (Pix2) (Proteome of centriole protein 1A) (WD repeat-containing protein 51A), POC1A, WDR51A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	POC1A

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA822753ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA822753ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.