





ZNF195 Antibody

Product Code	CSB-PA026570ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14628
Immunogen	Recombinant Human Zinc finger protein 195 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 195, ZNF195, ZNFP104
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF195

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA026570ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA026570ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.