





ZNF20 Antibody

Product Code	CSB-PA026573ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17024
Immunogen	Recombinant Human Zinc finger protein 20 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 20 (Zinc finger protein KOX13), ZNF20, KOX13
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF20

Image



Western blot
All lanes: ZNF20 antibody at 2.53μg/ml + MCF-7 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 62 kDa
Observed band size: 62, 40 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA026573ESR1HU at dilution of 1:100