




ZNF253 Antibody

Product Code	CSB-PA026613ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75346
Immunogen	Recombinant Human Zinc finger protein 253 protein (80-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 253 (Bone marrow zinc finger 1) (BMZF-1) (Zinc finger protein 411), ZNF253, BMZF1 ZNF411
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF253

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA026613ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.