



CACNG4 Antibody

Product Code	CSB-PA883364ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBN1
Immunogen	Recombinant Human Voltage-dependent calcium channel gamma-4 subunit protein (208-327AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-dependent calcium channel gamma-4 subunit (Neuronal voltage-gated calcium channel gamma-4 subunit) (Transmembrane AMPAR regulatory protein gamma-4) (TARP gamma-4), CACNG4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CACNG4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.