





CALCOCO1 Antibody

Product Code	CSB-PA878956ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P1Z2
Immunogen	Recombinant Human Calcium-binding and coiled-coil domain-containing protein 1 protein (462-691AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding and coiled-coil domain-containing protein 1 (Calphoglin) (Coiled-coil coactivator protein) (Sarcoma antigen NY-SAR-3), CALCOCO1, KIAA1536
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CALCOCO1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA878956ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer tissue using CSB-PA878956ESR2HU at dilution of 1:100