




CAMK1D Antibody

Product Code	CSB-PA812863ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IU85
Immunogen	Recombinant Human Calcium/calmodulin-dependent protein kinase type 1D protein (1-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium/calmodulin-dependent protein kinase type 1D (EC 2.7.11.17) (CaM kinase I delta) (CaM kinase ID) (CaM-KI delta) (CaMKI delta) (CaMKID) (CaMKI-like protein kinase) (CKLIk), CAMK1D, CAMKID
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CAMK1D

Image



Western blot

All lanes: CAMK1D antibody at 2.9µg/ml

Lane 1: A549 whole cell lysate

Lane 2: 293T whole cell lysate

Lane 3: Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 43, 41 kDa

Observed band size: 43 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.