





CEP120 Antibody

Product Code	CSB-PA836696ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N960
Immunogen	Recombinant Human Centrosomal protein of 120 kDa protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Centrosomal protein of 120 kDa (Cep120) (Coiled-coil domain-containing protein 100), CEP120, CCDC100
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CEP120

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA836696ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA836696ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.