





CXXC4 Antibody

Product Code	CSB-PA875656ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2H0
Immunogen	Recombinant Human CXXC-type zinc finger protein 4 protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	CXXC-type zinc finger protein 4 (Inhibition of the Dvl and axin complex protein), CXXC4, IDAX
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CXXC4

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA875656ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA875656ESR2HU at dilution of 1:100