




# DHDDS Antibody

<b>Product Code</b>	CSB-PA006845ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q86SQ9
<b>Immunogen</b>	Recombinant Human Dehydrodolichyl diphosphate synthase complex subunit DHDDS protein (1-333AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Dehydrodolichyl diphosphate synthase complex subunit DHDDS (EC 2.5.1.87) (Cis-isoprenyltransferase) (CIT) (Cis-IPTase) (Epididymis tissue protein Li 189m), DHDDS, HDS
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	DHDDS

## Image



Western blot


All lanes: DHDDSantibody at 2.97µg/ml + HeLa whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 39, 35 kDa

Observed band size: 39 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA006845ESR1HU at dilution of 1:100

---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.