




DHDDS Antibody

Product Code	CSB-PA006845ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86SQ9
Immunogen	Recombinant Human Dehydrodolichyl diphosphate synthase complex subunit DHDDS protein (1-333AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dehydrodolichyl diphosphate synthase complex subunit DHDDS (EC 2.5.1.87) (Cis-isoprenyltransferase) (CIT) (Cis-IPTase) (Epididymis tissue protein Li 189m), DHDDS, HDS
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DHDDS

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA006845ESR2HU at dilution of 1:100