





FGD5 Antibody

Product Code	CSB-PA757939ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZNL6
Immunogen	Recombinant Human FYVE, RhoGEF and PH domain-containing protein 5 protein (1203-1462AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	FYVE, RhoGEF and PH domain-containing protein 5 (Zinc finger FYVE domain-containing protein 23), FGD5, ZFYVE23
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FGD5

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA757939ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human bladder cancer using CSB-PA757939ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.