




GLYCTK Antibody

Product Code	CSB-PA009563ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVS8
Immunogen	Recombinant Human Glycerate kinase protein (284-523AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycerate kinase (EC 2.7.1.31) (HBeAg-binding protein 4), GLYCTK, HBEBP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GLYCTK

Image



Western blot

All lanes: GLYCTK antibody at 5.43µg/ml

Lane 1: Mouse liver tissue

Lane 2: Mouse lung tissue

Lane 3: Mouse kidney tissue

Lane 4: Rat lung tissue


Lane 5: NIH/3T3 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 56, 26, 49, 39, 20, 23 kDa


Observed band size: 56 kDa



Immunohistochemistry of paraffin-embedded

human lung tissue using CSB-

PA009563ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA009563ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.