





INHBC Antibody

Product Code	CSB-PA011721ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55103
Immunogen	Recombinant Human Inhibin beta C chain protein (103-352AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inhibin beta C chain (Activin beta-C chain), INHBC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INHBC

Image



Western blot
 All lanes: INHBC antibody at 5.16µg/ml
 Lane 1: Hela whole cell lysate
 Lane 2: A549 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 38 kDa
 Observed band size: 38 kDa



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA011721ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.