



ITPKB Antibody

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA011910ESR1HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P27987 |
| Immunogen | Recombinant Human Inositol-trisphosphate 3-kinase B protein (1-260AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Inositol-trisphosphate 3-kinase B (EC 2.7.1.127) (Inositol 1,4,5-trisphosphate 3-kinase B) (IP3 3-kinase B) (IP3K B) (InsP 3-kinase B), ITPKB |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | ITPKB |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |