





PIP5K1B Antibody

Product Code	CSB-PA018026ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14986
Immunogen	Recombinant Human Phosphatidylinositol 4-phosphate 5-kinase type-1 beta protein (190-470AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol 4-phosphate 5-kinase type-1 beta (PIP5K1-beta) (PtdIns(4)P-5-kinase 1 beta) (EC 2.7.1.68) (Phosphatidylinositol 4-phosphate 5-kinase type I beta) (PIP5KIbeta) (Protein STM-7) (Type I phosphatidylinositol 4-phosphate 5-kinase beta), PIP5K1B, STM7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIP5K1B

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA018026ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA018026ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.