





RAB17 Antibody

Product Code	CSB-PA863944ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0T7
Immunogen	Recombinant Human Ras-related protein Rab-17 protein (1-212AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-17, RAB17
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAB17


Image



Western blot
All lanes: RAB17 antibody at 3.02μg/ml
Lane 1: Mouse brain tissue
Lane 2: Mouse kidney tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 24, 10 kDa
Observed band size: 24 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA863944ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA863944ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.