




TM9SF3 Antibody

Product Code	CSB-PA023628ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HD45
Immunogen	Recombinant Human Transmembrane 9 superfamily member 3 protein (133-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane 9 superfamily member 3 (EP70-P-iso) (SM-11044-binding protein), TM9SF3, SMBP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TM9SF3

Image


Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA023628ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.