



TNNI3K Antibody

Product Code	CSB-PA701578ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q59H18
Immunogen	Recombinant Human Serine/threonine-protein kinase TNNI3K protein (1-75AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase TNNI3K (EC 2.7.11.1) (Cardiac ankyrin repeat kinase) (Cardiac troponin I-interacting kinase) (TNNI3-interacting kinase), TNNI3K, CARK
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TNNI3K

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA701578ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.