




TNNI3K Antibody

Product Code	CSB-PA701578ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q59H18
Immunogen	Recombinant Human Serine/threonine-protein kinase TNNI3K protein (1-75AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase TNNI3K (EC 2.7.11.1) (Cardiac ankyrin repeat kinase) (Cardiac troponin I-interacting kinase) (TNNI3-interacting kinase), TNNI3K, CARK
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TNNI3K

Image



Western blot

All lanes: TNNI3K antibody at 4.77μg/ml

Lane 1: Hela whole cell lysate

Lane 2: Mouse heart tissue

Lane 3: HepG2 whole cell lysate


Lane 4: Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 93, 105, 78, 95 kDa

Observed band size: 93 kDa



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA701578ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.