





UGT3A2 Antibody

Product Code	CSB-PA667485ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3SY77
Immunogen	Recombinant Human UDP-glucuronosyltransferase 3A2 protein (23-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-glucuronosyltransferase 3A2 (UDPGT 3A2) (EC 2.4.1.17), UGT3A2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UGT3A2

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA667485ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA667485ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.