




USP26 Antibody

Product Code	CSB-PA874835ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXU7
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 26 protein (664-913AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 26 (EC 3.4.19.12) (Deubiquitinating enzyme 26) (Ubiquitin thioesterase 26) (Ubiquitin-specific-processing protease 26), USP26
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP26

Image



Western blot

All lanes: USP26 antibody at 5.22μg/ml

Lane 1: Hela whole cell lysate

Lane 2: U937 whole cell lysate

Lane 3: Mouse gonadal tissue


Lane 4: Mouse liver tissue

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 105 kDa

Observed band size: 105 kDa



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA874835ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA874835ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.