




VANGL2 Antibody

Product Code	CSB-PA890750ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULK5
Immunogen	Recombinant Human Vang-like protein 2 protein (1-110AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vang-like protein 2 (Loop-tail protein 1 homolog) (Strabismus 1) (Van Gogh-like protein 2), VANGL2, KIAA1215 STB1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	VANGL2

Image



Western blot

All lanes: VANGL2 antibody at 2.26μg/ml

Lane 1: HepG2 whole cell lysate

Lane 2: MCF-7 whole cell lysate

Lane 3: Hela whole cell lysate


Lane 4: NIH/3T3 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 60 kDa


Observed band size: 60 kDa



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA890750ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA890750ESR1HU at dilution of 1:100



Western Blot

Positive WB detected in: HepG2 whole cell lysate, MCF-7 whole cell lysate, Mouse kidney tissue

All lanes: Vangl2 antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 60 kDa

Observed band size: 60 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.