





VANGL2 Antibody

Product Code	CSB-PA890750ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULK5
Immunogen	Recombinant Human Vang-like protein 2 protein (1-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vang-like protein 2 (Loop-tail protein 1 homolog) (Strabismus 1) (Van Gogh-like protein 2), VANGL2, KIAA1215 STB1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	VANGL2

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA890750ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA890750ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.