





WNT5B Antibody

Product Code	CSB-PA872432ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1J7
Immunogen	Recombinant Human Protein Wnt-5b protein (120-359AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Wnt-5b, WNT5B
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	WNT5B

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA872432ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA872432ESR1HU at dilution of 1:100