





DAZ2 Antibody

Product Code	CSB-PA623785ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13117
Immunogen	Recombinant Human Deleted in azoospermia protein 2 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deleted in azoospermia protein 2, DAZ2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DAZ2

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA623785ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA623785ESR1HU at dilution of 1:100