





SLC22A11 Antibody

Product Code	CSB-PA882097ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSA0
Immunogen	Recombinant Human Solute carrier family 22 member 11 protein (40-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 22 member 11 (Organic anion transporter 4), SLC22A11, OAT4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC22A11

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA882097ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA882097ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.