




ARL6 Antibody

Product Code	CSB-PA887945ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0F7
Immunogen	Recombinant Human ADP-ribosylation factor-like protein 6 protein (1-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylation factor-like protein 6 (Bardet-Biedl syndrome 3 protein), ARL6, BBS3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARL6

Image



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA887945ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.