




CCP110 Antibody

Product Code	CSB-PA528588ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43303
Immunogen	Recombinant Human Centriolar coiled-coil protein of 110 kDa protein (1-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Centriolar coiled-coil protein of 110 kDa (Centrosomal protein of 110 kDa) (CP110) (Cep110), CCP110, CEP110 CP110 KIAA0419
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCP110

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA528588ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.