




BVES Antibody

Product Code	CSB-PA002887ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NE79
Immunogen	Recombinant Human Blood vessel epicardial substance protein (241-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Blood vessel epicardial substance (hBVES) (Popeye domain-containing protein 1) (Popeye protein 1), BVES, POP1 POPDC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BVES

Image



Immunofluorescent analysis of 293 cells using CSB-PA002887ESR1HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.