





NUAK2 Antibody

Product Code	CSB-PA880931ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H093
Immunogen	Recombinant Human NUAK family SNF1-like kinase 2 protein (479-628AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NUAK family SNF1-like kinase 2 (EC 2.7.11.1) (Omphalocele kinase 2) (SNF1/AMP kinase-related kinase) (SNARK), NUAK2, OMPHK2 SNARK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NUAK2

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA880931ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA880931ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.