





NEK8 Antibody

Product Code	CSB-PA801219ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86SG6
Immunogen	Recombinant Human Serine/threonine-protein kinase Nek8 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase Nek8 (EC 2.7.11.1) (Never in mitosis A-related kinase 8) (NimA-related protein kinase 8) (Nima-related protein kinase 12a), NEK8, JCK NEK12A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NEK8

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA801219ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA801219ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.