



CST8 Antibody

Product Code	CSB-PA006096ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60676
Immunogen	Recombinant Human Cystatin-8 protein (53-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cystatin-8 (Cystatin-related epididymal spermatogenic protein), CST8, CRES
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CST8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.