





GNAO1 Antibody

Product Code	CSB-PA009593DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09471
Immunogen	Recombinant Human Guanine nucleotide-binding protein G(o) subunit alpha protein (60-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein G(o) subunit alpha, GNAO1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GNAO1

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA009593DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA009593DSR2HU at dilution of 1:100