




EEF1B2 Antibody

Product Code	CSB-PA335050DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24534
Immunogen	Recombinant Human Elongation factor 1-beta protein (1-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongation factor 1-beta (EF-1-beta), EEF1B2, EEF1B EF1B
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EEF1B2

Image



Western blot

All lanes: EEF1B2 antibody at 0.91 µg/ml

Lane 1: PC-3 whole cell lysate

Lane 2: 293T whole cell lysate

Lane 3: Hela whole cell lysate

Lane 4: Jurkat whole cell lysate


Lane 5: Raw264.7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 25 kDa

Observed band size: 25 kDa



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA335050DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.