




SLC1A5 Antibody

Product Code	CSB-PA620892DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15758
Immunogen	Recombinant Human Neutral amino acid transporter B(0) protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutral amino acid transporter B(0) (ATB(0)) (Baboon M7 virus receptor) (RD114/simian type D retrovirus receptor) (Sodium-dependent neutral amino acid transporter type 2) (Solute carrier family 1 member 5), SLC1A5, ASCT2 M7V1 RDR RDRC
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	SLC1A5

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA620892DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.