




# BCS1L Antibody

|                            |                                                                        |
|----------------------------|------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA897463ESR2HU                                                     |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.          |
| <b>Uniprot No.</b>         | Q9Y276                                                                 |
| <b>Immunogen</b>           | Recombinant Human Mitochondrial chaperone BCS1 protein (33-240AA)      |
| <b>Raised In</b>           | Rabbit                                                                 |
| <b>Species Reactivity</b>  | Human, Mouse                                                           |
| <b>Tested Applications</b> | ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200 |
| <b>Form</b>                | Liquid                                                                 |
| <b>Conjugate</b>           | Non-conjugated                                                         |
| <b>Storage Buffer</b>      | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                      |
| <b>Purification Method</b> | Antigen Affinity Purified                                              |
| <b>Isotype</b>             | IgG                                                                    |
| <b>Clonality</b>           | Polyclonal                                                             |
| <b>Alias</b>               | Mitochondrial chaperone BCS1 (h-BCS1) (BCS1-like protein), BCS1L, BCS1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                   |
| <b>Research Area</b>       | Signal Transduction                                                    |
| <b>Target Names</b>        | BCS1L                                                                  |

## Image



### Western blot

All lanes: BCS1L antibody at 5.04 µg/ml

Lane 1: Mouse kidney tissue

Lane 2: Mouse liver tissue

Lane 3: Mouse brain tissue

Lane 4: A549 whole cell lysate

Lane 5: 293T whole cell lysate


Lane 6: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48 kDa

Observed band size: 48, 30 kDa



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA897463ESR2HU at dilution of 1:100

---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.