




ERLIN2 Antibody

Product Code	CSB-PA007791ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94905
Immunogen	Recombinant Human Erlin-2 protein (212-339AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Erlin-2 (Endoplasmic reticulum lipid raft-associated protein 2) (Stomatin-prohibitin-flotillin-HflC/K domain-containing protein 2) (SPFH domain-containing protein 2), ERLIN2, C8orf2 SPFH2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ERLIN2

Image



Western blot

All lanes: ERLIN2 antibody at 3.1µg/ml

Lane 1: 293T whole cell lysate

Lane 2: HepG2 whole cell lysate

Lane 3: U251 whole cell lysate


Lane 4: Hela whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 38, 17, 23 kDa

Observed band size: 38 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA007791ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA007791ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.