




FBXO7 Antibody

Product Code	CSB-PA667562ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3I1
Immunogen	Recombinant Human F-box only protein 7 protein (353-522AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 7, FBXO7, FBX7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXO7

Image



Western blot

All lanes: FBXO7 antibody at 3.97µg/ml

Lane 1: Mouse gonadal tissue

Lane 2: Jurkat whole cell lysate

Lane 3: 293T whole cell lysate


Lane 4: PC-3 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 59, 50, 46 kDa

Observed band size: 59 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA667562ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA667562ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.