




FMO5 Antibody

Product Code	CSB-PA008751ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49326
Immunogen	Recombinant Human Dimethylaniline monooxygenase [N-oxide-forming] 5 protein (1-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dimethylaniline monooxygenase [N-oxide-forming] 5 (EC 1.14.13.8) (Dimethylaniline oxidase 5) (Hepatic flavin-containing monooxygenase 5) (FMO 5), FMO5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FMO5

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA008751ESR1HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.